

ResMed's Non-vented (NV) Mask Series comprises full face and nasal masks, to cater for a range of users and therapy needs.

The Series—which includes Ultra Mirage NV nasal mask, and Quattro FX NV and Ultra Mirage NV full face masks—provides a secure and comfortable solution for noninvasive pressure therapy in the hospital or the home.

NV Nasal Masks

Ideal where mouth leak is not a significant concern for patients. The masks may also be a valuable part of treatment in the home, where patient comfort is important.



Ultra Mirage™ NV
NASAL MASK

NV Full Face Masks

Covering both the nose and mouth, these masks are an effective solution for mouth leak and maintaining ventilation, especially at higher pressure levels. Suitable for long-term noninvasive ventilation in the home as well as the hospital environment.



Quattro FX NV
FULL FACE MASK



Ultra Mirage™ NV
FULL FACE MASK

Look for the blue tint on the mask frame or elbow: a safety feature that easily distinguishes non-vented masks from vented ones.

Non-vented (NV)
MASK SERIES

Product Information



Quattro FX NV
Full Face Mask



Ultra Mirage NV
Full Face Mask



Ultra Mirage NV
Nasal Mask

Features	Quattro FX NV Full Face Mask	Ultra Mirage NV Full Face Mask	Ultra Mirage NV Nasal Mask
Mirage dual-wall cushion technology	■	■	■
Spring Air cushion and Spring frame	■	■	■
Clear (non-tinted) frame	■	■	■
Breath-O-Prene headgear	■	■	■
Mask frame sizes	S, M, L	S, M, L	One size fits most
Cushion sizes	S, M, L	Standard, Shallow	Standard, Large, Shallow, Shallow Wide
Total number of mask size options	3	6	4
Number of points connecting frame to headgear	4	4	4
Quick-release headgear clips	2	4	2
Emergency release cord		■	
Adjustable forehead support	N/A	■	
Flexible forehead pads	N/A	■	
Suitable for patients with mouth leak	■	■	
360° rotating elbow	■	■	■
Cushion clips		■	■
Fitting/sizing tools	Custom-made sizing tool and Swing Seal fitting technique	3D fitting template	Fitting template

Product Codes

Quattro FX NV Full Face Mask

Asia Pacific, ANZ	Europe 1	Europe 2	Europe 3	USA
Small 61755	Small 61742	Small 61749	Small 61752	Small 61739
Medium 61756	Medium 61743	Medium 61750	Medium 61753	Medium 61740
Large 61757	Large 61744	Large 61751	Large 61754	Large 61741

Ultra Mirage NV Full Face Mask

Europe 1, ANZ, Asia Pacific	Europe 2	Europe 3	USA	3D Fitting Template for Ultra Mirage NV FFM
Small Standard 60643	Small Standard 60648	Small Standard 60606	Small Shallow 60636	English 16981
Small Shallow 60642	Small Shallow 60649	Small Shallow 60607	Small Standard 60637	French 16982
Medium Shallow 60644	Medium Standard 60650	Medium Standard 60608	Medium Shallow 60638	German 16983
Medium Standard 60645	Medium Shallow 60651	Medium Shallow 60609	Medium Standard 60639	
Large Shallow 60646	Large Standard 60652	Large Standard 60616	Large Shallow 60640	
Large Standard 60647	Large Shallow 60653	Large Shallow 60617	Large Standard 60641	

Ultra Mirage NV Nasal Mask

Europe* and Asia Pacific	USA
Standard 60401	Standard 60404
Large 60402	Large 60405
Shallow 60403	Shallow 60406
Shallow Wide 50408	

Europe 1: English, German, French, Spanish, Italian, Portuguese and Dutch
 Europe 2: English, Swedish, Norwegian, Danish and Finnish
 Europe 3: English, German, French, Spanish, Italian, Portuguese, Dutch, Hungarian, Russian, Polish, Turkish, Czech, Slovakian and Greek
 * English (UK), French, German and Italian (Switzerland)

ResMed's Non-vented Mask Series



Quattro FX NV
FULL FACE MASK

Non-vented (NV)
MASK SERIES

Security and comfort for noninvasive ventilation

All ResMed non-vented masks carry ResMed's renowned quality, comfort and performance.



The Series

ResMed's non-vented masks are suitable for patients suffering from acute as well as chronic respiratory distress or failure. The masks are specifically designed to be connected to ventilators that use an active exhaust valve and/or breathing circuit with an exhaust valve.

As well as housing many innovative features, the masks are made using hypoallergenic (latex-free), durable materials that add to patient comfort and extend mask life; are easy to assemble, set up and adjust; and have low dead space, contributing to reduced CO₂ rebreathing.

Synchronization is key

Patient-ventilator synchronization is the key to effective noninvasive ventilation. Poor synchronization can be caused by mask leak, which affects the ability of ventilators to accurately detect the patient's breathing pattern. ResMed's renowned cushion technology optimizes seal and minimizes leak, resulting in ideal synchronization.

Quattro™ FX NV FULL FACE MASK



The unobtrusive solution

- **Clear line of sight**—a feeling of openness created by a full face mask that starts mid-way on the nasal bridge.
- **New technologies in the Spring Air™ cushion and Spring frame** offer the stability and angular adjustment provided by conventional forehead supports.
- **It's as easy as Size. Swing. Seal**—Quattro FX NV gives you the tools to achieve the best fit for your patients. A custom-made sizing tool and Swing Seal fitting technique simplify the sizing, fitting and sealing process.
- **Contoured cushion membrane around sensitive nasal bridge area**—rolls and inflates to gently meet the nose rather than stretch over it.
- **New, tension-reducing headgear design**—features a contoured crown strap that sits snug at the top of the head and avoids the sensitive area around the back of the neck.

Ultra Mirage™ NV FULL FACE MASK



The conventional solution

- **Polycarbonate Luer lock port caps**—easier to use than conventional silicone caps.
- **360° rotating elbow**—provides flexibility and control over the breathing circuit, and optimises usability, comfort and mobility.
- **Unique adjustable forehead support, with four positions**—shortens fitting time for each individual patient and helps maintain seal even during slight patient movements, such as opening the mouth. The most effective position can be selected while the patient is wearing the mask.
- **Flexible forehead pads**—conform to the contour of the patient's forehead. The pads distribute pressure evenly over a larger surface area, reducing pressure on the face.
- **Quick-release headgear clips**—allow mask to be removed easily, without losing optimal headgear settings.
- **Emergency release cord**—allows patients to remove the mask in one simple action in an emergency.

Ultra Mirage™ NV NASAL MASK



The nasal solution

- **360° rotating elbow**—provides control over tubing position and breathing circuit, and optimizes usability, comfort and mobility.
- **Three cushion sizes**—large, standard and shallow provide more fitting options to suit more patients.
- **Forehead support**—easy, one-handed adjustment with four positions to maximize seal and patient comfort.
- **Four points of attachment between mask frame and headgear**—enhance stability by reducing possible mask movement.
- **Cushion clip**—securely fastens the cushion to the frame, improving strength, security and durability.
- **Quick-release headgear clips**—easy detachment in an emergency and easy fitting, even when the patient is in bed; also minimize the need for continual adjustment of the tabs.

Masks can be cleaned and reused between patients—all mask components can be disinfected/sterilized using validated cleaning procedures. Refer to the ResMed Masks Disinfection/Sterilization Guide online for details. www.resmed.com/masks/sterilization

Non-vented (NV)
MASK SERIES